

# Chapter 16

## Imaging of Cellular Oxygen and Analysis of Metabolic Responses of Mammalian Cells

Andreas Fercher, Tomas C. O’Riordan, Alexander V. Zhdanov,  
Ruslan I. Dmitriev, and Dmitri B. Papkovsky

### Abstract

Many parameters reflecting mitochondrial function and metabolic status of the cell, including the mitochondrial membrane potential, reactive oxygen species, ATP, NADH, ion gradients, and ion fluxes ( $\text{Ca}^{2+}$ ,  $\text{H}^+$ ), are amenable for analysis by live cell imaging and are widely used in many labs. However, one key metabolite – cellular oxygen – is currently not analyzed routinely. Here we present several imaging techniques that use the phosphorescent oxygen-sensitive probes loaded intracellularly and which allow real-time monitoring of  $\text{O}_2$  in live respiring cells and metabolic responses to cell stimulation. The techniques include conventional wide-field fluorescence microscopy to monitor relative changes in cell respiration, microsecond FLIM format which provides quantitative readout of  $\text{O}_2$  concentration within/near the cells, and live cell array devices for the monitoring of metabolic responses of individual suspension cells. Step by step procedures of typical experiments for each of these applications and troubleshooting guide are given.

**Key words:** Cell metabolism, mitochondrial function, cellular oxygen, phosphorescent oxygen-sensitive probe, phosphorescence quenching, live cell imaging, microsecond FLIM.

---

### 1. Introduction

Molecular oxygen ( $\text{O}_2$ ) is the key substrate of aerobic organisms and the terminal acceptor of the electron transport chain, which can serve as an informative marker of cell metabolism and mitochondrial function (1–6). The development of a number of fluorescence and phosphorescence based  $\text{O}_2$ -sensitive probes and assays (7–10) has provided relatively simple, non-chemical and

non-invasive means to measure  $O_2$  in samples containing respiring cells. Based on the collisional (i.e., non-chemical) quenching by  $O_2$  of certain photoluminescent structures (8), these probes combine good selectivity and reliable work across the whole physiological range of  $O_2$  (0–250  $\mu M$ ) and particularly at low  $O_2$  concentrations. They allow contact-less “sensing” of  $O_2$  in complex biological samples and particularly in samples containing respiring mammalian cells. Using extracellular  $O_2$ -sensitive probes, optical  $O_2$  sensing and respirometry were initially applied to the imaging of tissue oxygenation (9) and analysis of bulk consumption of  $O_2$  by populations of cells. A panel of high-throughput  $O_2$  consumption assays performed in standard microtiter plates with detection on a fluorescent plate reader have been developed and successfully used in a number of useful applications (10, 11–15).

More recently,  $O_2$  probes suitable for intracellular use have been developed (16, 17–22), which enable the analysis of localized  $O_2$  gradients within the cells and a number of additional measurement tasks. In particular, real-time monitoring of fast, transient respiratory responses of individual cells and cell populations (16, 20, 23, 24), which thus far was outside the scope of traditional  $O_2$ -sensing probes and techniques, has been achieved. Some of these probes are compatible for fluorescence microscopy imaging. It has been shown that using these intracellular  $O_2$  probes and relatively simple cell-handling procedures, sensing of intracellular  $O_2$ , and real-time monitoring of metabolic and respiratory responses of live mammalian cells and tissue can be performed (15, 21, 23). Several different cell types, particularly PC12, HepG2, HCT116, HeLa, SH-SY5Y, and A549 cells were examined so far by  $O_2$  imaging, which has been shown to provide information-rich data about cell metabolism and good spatial and temporal resolution (16, 22).

Application of the  $O_2$  imaging approach to suspension cells has been limited due to cell mobility and difficulties with sample manipulation and effector addition. These limitations have been overcome by the use of Live Cell Array<sup>TM</sup> (LCA) – a convenient platform developed for imaging of suspension cells (24). LCA provides capturing of cells within an array of specially designed microwells, whereby each well retains a single cell while maintaining cell viability. LCAs facilitate high content image analysis of large number of individual suspension cells with effector treatments.

In this chapter we describe the use of intracellular  $O_2$ -sensing probes based on phosphorescent platinum(II)-porphyrin dyes, in conjunction with fluorescence imaging for the assessment of mitochondrial function and real-time monitoring of responses of mammalian cells to metabolic effectors. These probes have been developed in our lab and commercialized by Luxcel Biosciences (Cork, Ireland). Standard wide-field fluorescence intensity

imaging in time-lapse mode, as well as microsecond FLIM with pulsed LED excitation are described and applied to commonly used cell models including the adherent HCT116 cells, differentiated PC12 cells, and suspension Jurkat T cells.

---

## 2. Materials

1. HCT116 human colorectal carcinoma, PC12 rat pheochromocytoma, and Jurkat human T-lymphoma cell lines (American Tissue Culture Collection, ATCC) (*see Note 1*).
2. Cell culture media: McCoy's 5A modified medium (McCoy), Roswell Park Memorial Institute 1640 (RPMI) (Sigma-Aldrich), Dulbecco's modified Eagle's medium (DMEM) (Invitrogen) (*see Note 2*).
3. Fetal bovine serum (FBS) and horse serum (HS) (Sigma).
4. Tissue culture grade dimethylsulfoxide (DMSO) and ethanol.
5. Penicillin/streptomycin (antibiotic/antimycotic) stock solution (Sigma)
6. Collagen IV and poly-D-lysine (Sigma) – prepared as 5  $\mu\text{g}/\text{mL}$  solution in growth medium for coating the dishes for cell adhesion.
7. Nerve growth factor (NGF) (Sigma) for differentiation of PC12 cells.
8. MitoXpress<sup>TM</sup> and IC60N phosphorescent oxygen-sensitive probes (Luxcel Biosciences, Ireland) (*see Note 3*).
9. EndoPorter<sup>TM</sup> reagent (Gene Tools, USA) for probe transfection – prepared as 1 mM stock in DMSO (*see Note 4*).
10. Glass-bottom imaging dishes (MatTek, USA).
11. Custom-made glass cylinders, approximately 9 mm outer diameter, 7 mm inner diameter, 15 mm long – one per dish (*see Note 5*).
12. Vacuum grease Apiezon L (Apiezon Corporation, USA).
13. Metabolic effectors (Sigma): mitochondrial uncoupler carbonylcyanide-*p*-trifluoromethoxyphenylhydrazone (FCCP) and complex III inhibitor antimycin A – prepared as 1 mM stock solutions in DMSO; extracellular  $\text{Ca}^{2+}$  chelator ethylene glycol-bis(2-aminoethylether)-N,N,N',N'-tetraacetic acid (EGTA) – prepared as 0.5 M solution in water, pH 7.4.

14. Live Cell Array<sup>TM</sup> devices (Nunc) for imaging of suspension cells.
15. Standard tissue culture equipment: CO<sub>2</sub> incubator, laminar flow cabinet, Gilson pipettes, centrifuge, vacuum aspirator, inverted transmission microscope, hemocytometer.

## **2.1. Imaging Equipment**

1. Standard wide-field fluorescent microscope Axiovert 200 with heated stage and objective (Carl Zeiss, Germany).
2. Fast gated CCD camera (LaVision, Germany) coupled to the microscope.
3. Custom-made LED excitation module (LaVision) coupled to the microscope, synchronized with CCD camera and controlled by the microscope software (LaVision). The module contains LED driver and a set of bright 397, 470, and 590 LEDs which are selected manually.
4. Filter cubes for the O<sub>2</sub>-probe: 395 nm bandpass excitation, 450 nm dichroic mirror, 650 nm longpass emission (Carl Zeiss) (*see Note 6*). Filter cubes for the other fluorescent probes (e.g., 488/525, 535/550 nm – Carl Zeiss), as required.
5. Incubation system for cell culture (PeCon, Germany) consisting of a heating insert mounted on the microscope stage, low-volume incubation chamber, temperature, and O<sub>2</sub> and CO<sub>2</sub> controllers. This system allows rapid adjustment of temperature, gas composition (O<sub>2</sub>, CO<sub>2</sub>), and humidity in the sample compartment and changing of O<sub>2</sub> concentration.
6. Perfusion pump (Ismatec SA, Glattbrugg, Switzerland) with 0.02 mm i.d. tubing (Gradko International Limited, Winchester, England) and a standard 22.5G hypodermic needle which are connected to the incubation chamber (*see details in Section 16.3*). The set-up provides effector addition during time-lapse imaging experiments, without disturbing the sample, gas composition, temperature, and focus of the region of interest (ROI).
7. ImSpector software (LaVision) for instrument control, FLIM and time-lapse measurements, image analysis, and data processing.
8. Wide-field fluorescence microscope IX51 (Olympus) equipped with a long distance 40× objective LUCPFLN (Olympus), climate control chamber for imaging in Live Cell Array devices. Equipped with a 395 nm excitation filter and a 600 nm cut-off emission filter.

### 3. Methods

#### 3.1. Culturing and Preparing the Cells for Imaging Experiments

1. All cell culture work is performed under a biological safety cabinet. HCT116 cells are cultured as adherent cells; PC12 cells are initially cultured in suspension and then adhered on pre-coated surfaces and differentiated (*see* **Notes 1** and **2**). Jurkat cells are cultured in suspension.
2. HCT116 cells are cultured in 75 cm<sup>2</sup> cell culture flasks with 0.2 µm vent cap in McCoy's 5A medium containing 2 mM L-glutamine, 100 U/mL penicillin, 100 µg/mL streptomycin (P/S), 10% FBS, maintaining them in CO<sub>2</sub> incubator at 37°C and 5% CO<sub>2</sub>. Cells are grown to a density of  $1\text{--}2 \times 10^6$  cells/mL, regularly split, and reseeded at  $1.5 \times 10^5$  cells/mL in 10 mL of medium. Cells are harvested from the flask using trypsin/EDTA, counted on a haemocytometer, and plated on imaging dishes by dispensing 0.2 mL of cell suspension ( $6 \times 10^4$  cells/mL) in the medium. Cells are then grown in minidishes to a confluence of approximately 90% (*see* **Note 7**) changing growth medium every 2 days.
3. PC12 cells are cultured in 75-cm<sup>2</sup> flasks in RPMI medium supplemented with 2 mM L-glutamine, 10% horse serum (HS), 5% fetal bovine serum (FBS), P/S at 5% CO<sub>2</sub> with regular passages in 30 mL at  $0.5\text{--}1 \times 10^5$  cells/mL. For seeding, cells are harvested, centrifuged at 200g for 5 min, separated from the supernatant, washed in PBS, centrifuged again, and resuspended in 1 mL of trypsin/EDTA solution. Cells are then incubated for 90 s at 37°C, diluted with 10 mL of RPMI supplemented with 1% HS, and gently passed ten times through 23G needle to separate individual cells. The cells are counted and seeded in MatTek dishes pre-coated with a mixture of collagen IV (0.007%) and poly-D-lysine (0.003%) at  $2.5 \times 10^4$  cells/dish. When the cells reach the confluence of ~50%, they are differentiated for 3–5 days in RPMI supplemented with 1% horse serum, P/S, and 100 ng/mL NGF. After differentiation cell layer becomes confluent and cells remain active for 2–4 days (*see* **Notes 7** and **8**).
4. Jurkat cells are cultured in 75-cm<sup>2</sup> flasks in RPMI containing 10% FBS, P/S to a density of  $1\text{--}2 \times 10^6$  cells/mL at 37°C and 5% CO<sub>2</sub>. Cells are split regularly and reseeded at  $1.5 \times 10^5$  cells/mL in 20 mL of medium (*see* **Note 8**). Cells are then harvested from the flask, counted on a hemocytometer, centrifuged at 200g, resuspended in medium, adjusted to a concentration  $1 \times 10^6$  cells/mL, and kept at 37°C.

### **3.2. Imaging of Adherent HCT116 Cells Loaded with MitoXpress Probe**

1. Reconstitute a 1 × package of MitoXpress probe (dry sample in a plastic vial) in 0.1 mL of growth medium with additives used to culture HCT116 cells, to produce a 10- $\mu$ M probe stock (*see Note 9*).
2. Pipette the required amount of medium (0.2 mL per dish) to a glass or plastic vial, add stock solutions of probe (100  $\mu$ L/mL) and Endoport<sup>®</sup> (6  $\mu$ L/mL), and incubate the resulting probe loading solution for 10–15 min at 37°C.
3. Take the dishes with HCT116 cells prepared for imaging experiments as described in **Section 3.1**, step 2, and transfer 200  $\mu$ L of probe loading solution to each dish.
4. Incubate the dishes in CO<sub>2</sub> incubator at 37°C for 24–30 h (*see Note 10*).
5. After the incubation, carefully aspirate probe loading solution from the dish and discard it. Add 1 mL of fresh medium, incubate for 5 min, then aspirate and discard. Repeat washing of the cells two more times and finally add 0.2 mL of medium (*see Note 11*).
6. Switch on the microscope and cell incubation system and allow them to warm up and equilibrate (T, CO<sub>2</sub>, O<sub>2</sub>) for at least 30 min. Set up the software for image acquisition and time-lapse experiments and the required filter combination (395 nm excitation and >650 nm emission).
7. Take clean glass cylinder and apply a thin, uniform layer of vacuum grease to one of its edges.
8. Take the imaging dish with cells (prepared in **Section 3.1**), aspirate medium, and then firmly attach the cylinder with its greased side to the glass bottom of the dish where the cells are grown (*see Note 12*). Add 180  $\mu$ L of medium to the cap and 1 mL of medium to the outer space and place the measurement unit inside the heating insert. The assembled unit with glass cap is shown in **Fig. 16.1**.
9. Assemble the incubation chamber and allow the dish with the cap and cells to equilibrate for about 20 min (*see Note 13*).
10. Using a 100- $\mu$ L Hamilton syringe, load the solution of effector into the tubing and pump it inside the chamber to allow temperature and gas equilibration.
11. Switch the microscope in transmission mode; focus on cells inside the cap area.
12. Switch the microscope to fluorescence mode. Select ROI where the cells have relatively uniform density and good loading with probe. Representative bright field and




Fig. 16.1. Photographs of the MatTek dish with cap used for  $O_2$  imaging of adherent cells. (a) General view of 35-mm imaging dish with attached glass cylinder mounted using vacuum grease. (b) Incubation chamber of the microscope with a sample prepared for  $O_2$  measurements: 1 – glass cylinder; 2 – imaging dish; 3 – syringe needle connected to the pump for application of effector solution; 4 – heating insert with sample holder of the incubation system; 5 – microscope objective. Measurement requires restricted surface area to reduce convection and mass exchange within the sample, which is provided by a cylindrical cap attached to the bottom of the dish with cells. Effector addition should not disturb the sample. Small volume of effector stock solution is applied on top of the medium by means of perfusion pump with bended needle or gently applied onto the wall of the cap with a micropipette.

fluorescent images of HCT116 cells loaded with probe are shown in **Fig. 16.2a, b**.

13. Initiate image acquisition in time-lapse mode and monitor basal fluorescent signal of the cells under resting conditions.
14. If the signals are sufficiently high (compared to blanks outside the cells) and stable, proceed to the experiments with effector treatments (*see Note 14*). Initially, it is recommended to test cellular response by treating the cells with uncoupler FCCP and/or with mitochondrial inhibitor Antimycin A (*see Note 15*).
15. Move the microscope stage to a new field with cells within the dish, adjust the focus, and initiate image acquisition in time-lapse mode.
16. After recording the basal signal of resting cells for 5–15 min, make the addition of FCCP between two measurements. We normally add 1/10 of the total sample volume in the cap (i.e., 20  $\mu\text{L}$ –180  $\mu\text{L}$ ) using the lowest possible flow rate of the pump (*see Note 16*). Continue image acquisition for 10–30 min after the addition of FCCP. When the response is finished, stop the acquisition and analyze signal profile(s).
17. Interpretation of signal profiles of the intracellular probe. Basal fluorescent signal prior to effector addition reflects  $O_2$  concentration within the cell which corresponds to the resting level of respiration. An increase in probe signal after FCCP addition reflects increased respiration which




Fig. 16.2. Bright field (a) and fluorescent (b) images of HCT116 cells prepared for imaging in a MatTek dish and loaded with IC60 probe. (c) Respiratory response of HCT116 cells to the addition of FCCP measured by conventional imaging (390 nm bandpass/650 nm longpass). The increase in fluorescent signal reflects increased respiration due to mitochondrial uncoupling.

causes a reduction in cellular  $O_2$ . Conversely, Antimycin A is expected to inhibit cell respiration and elevate cellular  $O_2$  to the levels present in bulk medium (i.e., air-saturated), thus bringing probe signal to a lower level (*see Note 17*). A characteristic profile obtained with HCT116 cells and FCCP is shown in **Fig. 16.2c**.

18. If the responses to FCCP and Antimycin A treatments are clearly seen, the samples are prepared properly, and the cells are in good condition, one can proceed to the experiments with stimulants having unknown effects on cell respiration.
19. In this case, take a new dish with cells loaded with probe and follow steps 11–17 above. Make effector addition and



monitor relative changes in intracellular  $O_2$ , which reflect the changes in cellular respiration. If required, repeat treatment.

20. Quantitative assessment. Plot signal profile for ROI with cells; adjust intensity scale as required. Select ROI outside the cells and plot the profile of blank signal. Subtract blank profile from the cell signal profile to produce the profile of probe fluorescence (*see* **Note 18**).

### **3.3. FLIM Analysis of Responses of Differentiated PC12 Cells Loaded with IC60 Probe**

1. Prepare differentiated PC12 cells in imaging dishes as described in **Section 3.1**, step 3.
2. Take the vial containing IC60 probe, prepare stock solution and then load the cells with this probe using the same procedure as described above in **Section 3.2**, steps 1–5, for the MitoXpress probe.
3. Prepare samples and imaging system as described in **Section 3.2**, steps 6–10.
4. Start with a bright field image of the sample for pre-focusing, using repetition mode with an exposure time of 5 ms.
5. Change to fluorescence mode. Using steady-state mode (continuous excitation) and exposure time of 300–1000 ms, take one fluorescence image of the sample. Fluorescent signal received by the camera should be in the region of 1000–3000 counts per pixel (*see* **Note 19**). Make adjustments of the focus and sensitivity/scale. Repeat the adjustment taking more images if required.
6. Program the software for time-lapse experiments in the FLIM mode. Using ImSpector software, set up the parameters of phosphorescence lifetime measurements: excitation pulse width, camera modulation, gate time, and modulation frequency (ratio of the gate time and repetition time) (*see* **Note 20**). General scheme of phosphorescence decay measurement implemented in our system is shown in **Fig. 16.3**. Typical settings for the IC60 probe are given in **Table 16.1**.
7. Perform one set of measurements to generate phosphorescence decay curve for the ROI containing loaded cells. Make final adjustment of the sensitivity and intensity signals generated by FLIM. Intensity at the first delay time, which is a function of the exposure time, excitation pulse length, and x-y binning, should be in the range 1000–3000 counts (*see* **Note 21**). Signal intensity can be increased by binning; however, this will decrease image resolution.
8. Move the sample to a new field with cells and commence FLIM measurements. Each time measure




Fig. 16.3. Schematic description of FLIM measurements and the main parameters used: LED pulsing rate, gate time on the CCD camera, CCD camera exposure time, delay time, and the number of points on the curve (*left panel*). Measured phosphorescence decay curve (*right panel*).

**Table 16.1**  
**Typical settings for FLIM measurements with the O<sub>2</sub> probe**

Parameter	Setting
Exposure time	1–3 s
Pulse length	10 μs
Integration start time	2 μs after the pulse
Integration end time	Typically 60 μs
Number of measurement points	10–20 points
Gate time	10–30 μs
Cycle time	100 μs

phosphorescence decay of the whole image, repeating this every 2–3 min.

- Between two measurements of the phosphorescence decay make effector additions (refer to **Section 3.2**, steps 16 and 17 for details) and continue acquisition of decay curves. After the required number of measurements, stop the acquisition.
- Define ROIs within the sample which contain cells loaded with probe, and generate phosphorescence decay curves for each ROI and for each time point of kinetic measurement. Extrapolate these decay curves with single-exponential fits

and determine phosphorescence lifetime value for each measurement.

11. Plot lifetime values as a function of time for each of the ROI and analyze the resulting profiles. Examples of FLIM images and respiratory response of PC12 cells to the depletion of extracellular  $\text{Ca}^{2+}$  are shown in **Fig. 16.4**. Intensity profiles can also be generated from the images (*see Fig. 16.4b*). If required, phosphorescence lifetime profiles produced by FLIM can be converted into  $\text{O}_2$  concentration scale, using probe calibration which is determined separately (*see steps 12–16 below*).




Fig. 16.4. Respiratory response of differentiated PC12 cells loaded with IC60 probe to the depletion of extracellular  $\text{Ca}^{2+}$  measured by FLIM. (a) Distribution of the probe intensity signal and selected ROI (*inset*). (b) Profile of phosphorescence intensity within ROI showing signal-to-noise ratio. (c) The response to the addition of EGTA (2.5 mM) with a transient increase in probe lifetime reflecting partial deoxygenation of the cell due to increased respiration ( $p < 0.01$  at the peak of the response). (d) The relationship between probe lifetime and  $\text{O}_2$  concentration (calibration).

12. Take a fresh imaging dish (without glass cap) containing cells pre-loaded with the probe and 1–2 mL of medium and add Antimycin A to the medium at a final concentration of 10  $\mu\text{M}$ . This should block cellular respiration and eliminate  $\text{O}_2$  gradients in the sample.
13. Insert the dish in the incubation chamber with pre-set  $\text{pO}_2$  (e.g., 3%) and incubate it for about 30 min to achieve temperature and gas equilibrium.
14. Prepare the microscope for FLIM measurements as described above (*see steps 4–7 above*). When the sample is ready, measure phosphorescence decay curve. Calculate the lifetime as described in step 10 above.
15. Change the incubation system to a new  $\text{pO}_2$  setting and re-equilibrate the sample by incubating it for  $\sim 30$  min. Move the sample/microscope stage to a new field with cells and perform lifetime measurement again. Repeat this at 3–4 other  $\text{pO}_2$  settings.
16. Plot the relationship between phosphorescence lifetime of the probe loaded intracellularly and  $\text{O}_2$  concentration. Adjust it to the desired  $\text{O}_2$  scale: at  $37^\circ\text{C}$  and normal

ambient pressure, 20.5% O<sub>2</sub> in the chamber corresponds to 100% of air saturation or 200  $\mu$ M O<sub>2</sub>. Use the resulting relationship to convert measured phosphorescent lifetime profiles into profiles of cellular O<sub>2</sub> concentration.

### **3.4. Imaging in Live Cell Arrays of Suspension Jurkat Cells Loaded with the IC60 Probe**

1. Reconstitute 3  $\times$  package of IC60 probe (300 nmoles dried in a vial) in 0.3 mL of medium to produce a 10- $\mu$ M probe stock.
2. Take in a plastic vial 200  $\mu$ L of suspension of Jurkat cells ( $1 \times 10^6$  cells/mL) and add 20  $\mu$ L of probe and 1.2  $\mu$ L of Endoport stock solutions. Incubate at 37°C, 5% CO<sub>2</sub> for 24 h (*see Note 10*).
3. After incubation, wash the cells three times by centrifuging them at 200*g*, removing supernatant, and finally resuspend the cells in 0.1 mL of medium. Keep this suspension (approximately  $\sim 2 \times 10^6$  cells/mL) at 37°C, 5% CO<sub>2</sub> for further use.
4. Switch on the Olympus IX51 microscope. Take a new LCA device, place it on the stage of the Olympus microscope, and allow to warm up and equilibrate at 37°C for about 20 min (*see Note 13*).
5. Using a P20 pipette load the LCA with 2.5  $\mu$ L of suspension of cells prepared in step 3 above. Incubate for 5 min at 37°C and then flush with 15  $\mu$ L of fresh medium.
6. Using transmission mode, focus on the microwells with cells within the LCA. Ensure that the device is loaded properly and the majority of the microwells contain trapped cells (*see Note 22*). Sample images of the LCA with cells are shown in **Fig. 16.5a, b**.
7. Switch the microscope to fluorescence mode, adjust the focus and initiate image acquisition in time-lapse mode. Take images every 20–60 s for about 10 min (*see Note 14*), using excitation at 395 nm and collection of emission with a 600 nm cut-off filter. Make sure that cell loading with probe is sufficient, adjust the sensitivity and scale as required.
8. After basal respiration (resting cells) is recorded, pause monitoring and apply effector treatment. Take solution of effector in medium (diluted to the required final concentration) and apply with a Gilson P20 pipette 15  $\mu$ L of this solution to the loading bay of the LCA. After the effector is loaded in the LCA, resume monitoring for the required period of time (*see Note 23*).
9. Select several ROIs within the LCA – controls with cells – and analyze measured profiles and respiratory responses. If required, use the well without cells as a reference (blank).




Fig. 16.5. (a) Photograph of a working section of Live Cell Array unit used for imaging of suspension cells (reproduced from (24) with permission of the Royal Society of Chemistry). (b) Fluorescent image showing individual Jurkat cells loaded with the O<sub>2</sub> probe trapped in the cavities of LCA. (c) Metabolic response of individual Jurkat cell to the addition of mitochondrial uncoupler FCCP (1 μM), measured using the LCA.

For interpretation of measured signal profiles and respiratory responses of cells, refer to **Section 3.2**, step 18. Example of response of individual Jurkat cell to FCCP treatment is shown in **Fig. 16.5c**.

#### 4. Notes

1. A range of common cell lines can be used for imaging intracellular O<sub>2</sub>; however, probe loading efficiency and metabolic activity may vary greatly for different cell types. Cells having well-developed mitochondrial machinery and

high rates of O<sub>2</sub> consumption (e.g., primary cells) are the preferred choice.

2. Growth medium, additives, and optimal culturing conditions may vary for each particular cell type.
3. Other O<sub>2</sub> probes for intracellular use, including those which do not require transfection reagents, are also available or in development. Contact corresponding suppliers (e.g., Luxcel Biosciences, [www.luxcel.com](http://www.luxcel.com)) and research labs to get advice on the most suitable probe for your particular imaging system, measurement task, cell type.
4. In our experience, probe transfection with EndoPorter works well with many common cell types. Other loading reagents and procedures may be efficient as well. Generally, probe loading method and procedure need to be optimized for each cell type and medium.
5. Instead of the imaging dishes with cylinders, glass-bottom 96-well imaging microplates, multiwell flexiPERM inserts (Greiner Bio, Germany) can also be used. Devices with larger well size (e.g., 24-well plates) are less preferred.
6. Equivalent live cell fluorescent imaging systems, which are spectrally compatible and sensitive enough with the O<sub>2</sub> probes, can also be used. The MitoXpress<sup>TM</sup> probe is excitable at either 370–395 or 525–545 nm, and emits around 630–700 nm (maximum at 650 nm). The IC60 probe is excitable at 380–400 as well as at 580–600 nm, and it emits in the near infrared range 720–850 nm (maximum at 760 nm). IC60 probe is more photostable; however, its very longwave emission and lower brightness (compared to the MitoXpress<sup>TM</sup>) can be problematic for some imaging systems.
7. Relatively high surface density of cells and high respiratory activity are the main requirements of a successful metabolic assay using the intracellular O<sub>2</sub> probes. Cells having high density of mitochondria and high levels of respiration/oxidative phosphorylation such as PC12 and HCT116 produce robust responses to stimulation with various metabolic effectors, whereas slowly respiring and highly glycolytic cells would be difficult to analyze by this method.
8. Optimal conditions for culturing depend on particular cell type. Differentiation is important for many cells to display a range of functional characteristics observed in vivo.
9. Probe loading by facilitated endocytosis in the presence of Endoport reagent does not impose special requirements on the medium used. Other loading reagents and methods may require serum-free conditions.

10. Probe loading with EndoPorter is relatively slow; however, it works well in different media and cell types and produces high fluorescent signals which are easy to measure.
11. Extreme care should be taken not to damage the cell layer and wash the cells off the surface.
12. The cylinder is required to stabilize local oxygen gradients at cell layer by reducing surface area of the medium which in turn reduces convection and mass exchange within the sample. Based on our experience, this set-up produces better quality results.
13. Efficient temperature control of the measurement system and equilibration of the sample are very important. All O<sub>2</sub> probes produce fluorescence response to temperature fluctuations.
14. During these preparations, sample exposure to excitation light should be kept to a minimum. Photostability of the MitoXpress probe is moderate, whereas the IC60 probe is more photostable but less bright.
15. Optimal effector concentrations which produce maximal metabolic response depend on the cell type. We normally use 0.5–5  $\mu$ M for FCCP and 10  $\mu$ M for Antimycin A (in the well).
16. Mixing and disturbing the sample during the addition should be avoided as this will affect local O<sub>2</sub> gradients at cell monolayer and, hence, the probe signal.
17. Compared to mitochondrial inhibitors, uncouplers and activators of metabolism produce strong and positive optical response which is easier to measure. If basal respiration of cell monolayer is low and unable to create local O<sub>2</sub> gradient, metabolic responses may not be detectable, even though they are present. This is determined by the principle of operation of the intracellular O<sub>2</sub> assay (10). The lack or the absence of optical response from the sample/cells is a rather common case. This can be due to low respiration rate so that the cells are unable to create measurable local O<sub>2</sub> gradients and changes in response to stimulation. Cell damage during sample preparation, phototoxicity, or probe photobleaching during the measurement can also occur. The following can be used to improve assay sensitivity: (i) increase cell confluence (overgrow cells); (ii) increase cell respiration by using glucose-free medium containing 10 mM of galactose and pyruvate (12); (iii) increase the height of medium layer in the cap; (iv) conduct imaging experiments under reduced pO<sub>2</sub> (hypoxia chamber set at 10% or 5% pO<sub>2</sub>); (v) optimize cell culture and probe loading conditions.

18. If probe photobleaching is high, measured profiles should be corrected for it. The initial part of the profile, which corresponds to basal respiration, can be used to quantify the rate of photobleaching.
19. Try several different settings in steps 4 and 5, and select those which provide lowest photobleaching during image acquisition.
20. For our system, delay time is calculated from the ratio of the effective measurement time (between start and end of one cycle) and the number of measurement points less one point. Cycle time is the time between LED pulses. Gate is integration time of the fluorescence signal for each delay time.
21. To get approximately the same level of signal from the sample in the FLIM mode, exposure time on the camera needs to be increased approximately five times compared to the steady-state measurement mode.
22. If loading of the LCA with cells is poor, step 8 can be repeated.
23. The time of monitoring of metabolic responses of cells is determined by the mode of action of the effector. It is also limited by the changes in cell viability during the measurement (i.e., phototoxicity). Most of our imaging experiments were conducted over a period of 20–60 min. Since both MitoXpress and IC60 probes do not leak from the cells, longer experiments with reduced sampling frequency are possible.

---

## Acknowledgments

Financial support of this work by the European Commission (FP7 project CP-FP214706-2) and by the Science Foundation of Ireland grant (07/IN.1/B1804) is gratefully acknowledged.

## References

1. Aragonés J, Fraisl P, Baes M, Carmeliet P (2009). Oxygen sensors at the crossroad of metabolism. *Cell Metabol* **9**(1):11–22.
2. Duchen MR (2004). Mitochondria in health and disease: perspectives on a new mitochondrial biology. *Mol Aspects Med* **25**:365–451.
3. Dykens JA, Will Y (2007). The significance of mitochondrial toxicity testing in drug development. *Drug Discov Today* **12**:777–85.
4. Nicholls DG (2005). Mitochondria and calcium signaling. *Cell Calcium* **38**: 311–7.
5. Hynes J, Marroquin LD, Ogurtsov VI, Christiansen KN, Stevens GJ, Papkovsky DB, Will Y (2006). Investigation of drug-induced mitochondrial toxicity using fluorescence-based oxygen-sensitive probes. *Toxicol Sci* **92**:186–200.



6. Matoba S, Kang JG, Patino WD, Wragg A, Boehm M, Gavrilova O, Hurley PJ, Bunz F, Hwang PM (2006). p53 regulates mitochondrial respiration. *Science* **312**:1650–3.
7. Dunphy I, Vinogradov SA and Wilson DF (2002) Oxyphor R2 and G2: phosphors for measuring oxygen by oxygen-dependent quenching of phosphorescence. *Anal Biochem* **310**:191–198.
8. Papkovsky DB (2004). Methods in optical oxygen sensing: protocols and critical analyses. *Meth Enzymol* **381**:715–35.
9. Rumsey WL, Vanderkooi JM, Wilson DF (1988). Imaging of phosphorescence: a novel method for measuring oxygen distribution in perfused tissue. *Science* **241**:1649–51.
10. Hynes J, O’Riordan TC, Papkovsky DB (2008). The use of oxygen sensitive fluorescent probes for the assessment of mitochondrial function. In JA Dykens and Y Will (Eds.), *Mitochondrial Dysfunction in Drug Induced Toxicity*. Wiley & Son, New York.
11. Will Y, Hynes J, Ogurtsov VI, Papkovsky DB (2006). Analysis of mitochondrial function using phosphorescent oxygen-sensitive probes. *Nat Protoc* **1**(6):2563–72.
12. Marroquin LD, Hynes J, Dykens JA, Jamieson JD, Will Y (2007). Circumventing the Crabtree effect: replacing media glucose with galactose increases susceptibility of HepG2 cells to mitochondrial toxicants. *Toxicol Sci* **97**:539–47.
13. Wu M, Neilson A, Swift AL, Moran R, Tamagnine J, Parslow D, Armistead S, Lemire K, Orrell J, Teich J, Chomicz S, Ferrick DA (2007). Multiparameter metabolic analysis reveals a close link between attenuated mitochondrial bioenergetic function and enhanced glycolysis dependency in human tumor cells. *Am J Physiol Cell Physiol* **292**(1):C125–36.
14. Hynes J, O’Riordan TC, Zhdanov AV, Uray G, Will Y, Papkovsky DB (2009). In vitro analysis of cell metabolism using a long-decay pH-sensitive probe and extracellular acidification assay. *Anal Biochem* **390**(1):21–8.
15. O’Hagan KA, Cocchiglia S, Zhdanov AV, Tambawala MM, Cummins EP, Monfared M, Agbor TA, Garvey JF, Papkovsky DB, Taylor CT, Allan BB. (2009). PGC-1 $\alpha$  is coupled to HIF-1 $\alpha$ -dependent gene expression by increasing mitochondrial oxygen consumption in skeletal muscle cells. *Proc Natl Acad Sci USA* **106**(7):2188–93.
16. O’Riordan TC, Fitzgerald K, Ponomarev GV, Mackrill J, Hynes J, Taylor C, Papkovsky DB. (2007). Sensing intracellular oxygen using near-infrared phosphorescent probes and live-cell fluorescence imaging. *Am J Physiol Regul Integr Comp Physiol* **292**:R1613–20.
17. Ji J, Rosenzweig N, Jones I and Rosenzweig Z (2001). Molecular oxygen-sensitive fluorescent lipobeads for intracellular oxygen measurements in murine macrophages. *Anal Chem* **73**:3521–3527.
18. Xu H, Aylott JW, Kopelman R, Miller TJ and Philbert MA (2001). A real-time ratiometric method for the determination of molecular oxygen inside living cells using sol-gel-based spherical optical nanosensors with applications to rat C6 glioma. *Anal Chem* **73**:4124–4133.
19. Koo YE, Cao Y, Kopelman R, Koo SM, Brasuel M and Philbert MA (2004). Real-time measurements of dissolved oxygen inside live cells by organically modified silicate fluorescent nanosensors. *Anal Chem* **76**:2498–2505.
20. Schmalzlin E, van Dongen JT, Klimant I, Marmodee B, Steup M, Fisahn J, Geigenberger P, Lohmannsroben HG (2005). An optical multifrequency phase-modulation method using microbeads for measuring intracellular oxygen concentrations in plants. *Biophys J* **89**:1339–1345.
21. O’Riordan TC, Zhdanov AV, Ponomarev GV and Papkovsky DB (2007). Analysis of intracellular oxygen and metabolic responses of mammalian cells by time-resolved fluorometry. *Anal Chem* **79**:9414–9.
22. Finikova OS, Lebedev AY, Aprelev A, Troxler T, Gao F, Garnacho C, Muro S, Hochstrasser RM, Vinogradov SA (2008). Oxygen microscopy by two-photon-excited phosphorescence. *Chemphyschem* **9**(12):1673–9.
23. Zhdanov AV, Ward MW, Prehn JH and Papkovsky DB 2008. Dynamics of intracellular oxygen in PC12 Cells upon stimulation of neurotransmission. *J Biol Chem* **283**:5650–61.
24. Deutsch M, Deutsch A, Shirihai O, Hurevich I, Afrimzon E, Shafran Y, Zurgil N (2006). A novel miniature cell retainer for correlative high-content analysis of individual untethered non-adherent cells. *Lab Chip* **Aug**;6(8):995–1000.